







Features

- 1. High-quality poly surface for decubitus prevention, removable, washable, impermeable
- 2. Kidney bridge up
- 3. Backlit display indicating status information, e.g. currently chosen position, and status of battery charging, it also ensures a safe and easy handling even in the dark operation room.
- 4. Additional control at the bottom also is with backlit
- 5. Stainless steel inlay plates
- 6. Diam 125mm castors with locking, universal castors
- 7. Many selections of accessories, extensions, and positioning aids

Technical Specifications

Total length(with additional footplate): 2220mm

Total width: 600mm Net Width: 555mm

Length of Head section: 300mm

Length of Back section: 380mm

Length of Main section: 1080mm

Length of Leg section: 580mm

Supporting bottom: 910*625*140mm Height adjustment range: 520 -1000mm

Kidney support up: 120mm

Movement backward to forward: 350mm

Back section up: 70°







Technical Specifications

Back section down: 22° Trendelenburg: 35°

Anti-trendelenburg: -35° Lateral adjustment : 22°

Leg section inclination: +15° till -90°
Spread angle of leg plates: 100 degree
Head section inclination: +45° till -90°
Net weight of overall table: 280kgs
Patient horizontal max load: 300kgs

Max load weight: 460kgs

Integrated battery (back-up): For over 200 operations

Tables sections: 5 Radiolucent Top: Yes C-arm compatible: Yes

Foot: 40 degree Remote control: Yes Castor Lock: yes

Backlit display indicating status information, e.g. currently chosen position, and status of battery charging, it also ensures a safe and easy handling even in the dark operation room.



High quality poly surface for decubitus prevention, removable, washable, impermeable





Additional control at the bottom also is with backlit.



Algontech Canada Inc.

Email: info@algontech.com Phone: +1 6478755424 Website: algontech.com







Good performance



Stainless steel inlay plates





Diam 125mm castors with locking, universal castors

Diam 125mm castors, directional castors

Algontech Canada Inc.



Positions







Gynecology position

Nr.2 Gynecology / Giving birth

Endoscope application

Endoscope application

Nr.5 Neurosurgery

Neurosurgery in beach chair position

Neurosurgery













Nr.8 **ENT/ Neurosurgery**

Algontech Canada Inc.




Positions









Nr.10 Spinal surgery

Nr.11 Kidney position

Nr.12 General position (Table top movement

Nr.13 Urology/TUR

for C-arm)

Nr.14 Gynecology/Giving birth

Rectal surgery

Knee arthroscopy

















Positions







Nr.18 Upper limb trauma

Nr.19 Fractures of the neck of femur

Tibial shaft fracture








Algontech Canada Inc.





Basic features

















Algontech Canada Inc.





Basic features















Algontech Canada Inc.



Basic features



Algontech Canada Inc.

